

Smoke Evacuation Filters

Features:

- For passive removal of toxic smoke during laparoscopic surgery
- Filter connects to the gas outlet of the trocar
- Open the thumbwheel and smoke is filtered and released

Specifications:


- Viral/Bacterial Efficiency: 99.999%
- Deodorization Efficiency: 48%
- Minimum Particulate Size: 0.26 μ
- Filtration Surface Area: 33cm²
- Typical Air Flow Rate: 15L/min @ 1psi
- Maximum Operating Pressure: 30psi

Presentation:

- Sterile, single use

Sterilisation: Ethylene Oxide

ARTG No: 152979



PRODUCT CODE

AMD-ASF-3000	Smoke Filter, Sterile, Single Use, Box 25
--------------	---



For more information, call
Customer Service on 1300 616 366
or your local Product Specialist