

## Matrix Wound Retractors

### Features:

- Medical grade TPU
- Protects local tissue during surgery
- Providing a cleaner surgical environment
- Numerous sizes available
- Applications: Open Surgery - Abdominal, Colorectal, Gynae (Cesarean), ENT (Thyroidectomy)

### Presentation

- Each, in blister pouch

Sterilisation: Ethylene Oxide

ARTG No: 275852



### PRODUCT CODES

	Outer ring	Inner ring	Channel	Incision size
<b>IN2-WR-B1010</b>	60mm	70mm	60mm x 150mm	20-40mm
<b>IN2-WR-B1011</b>	80mm	90mm	80mm x 150mm	30-60mm
<b>IN2-WR-B1020</b>	120mm	130 mm	120mm x 250mm	50-100mm
<b>IN2-WR-B1021</b>	150mm	160mm	150mm x 250mm	80-130mm
<b>IN2-WR-B1022</b>	180mm	190mm	180mm x 250mm	100-150mm
<b>IN2-WR-C1023</b>	200mm	210mm	200mm x 250mm	80-150mm
<b>IN2-WR-C1030</b>	220mm	230mm	220mm x 250mm	80-150mm
<b>IN2-WR-C1031</b>	250mm	260mm	250mm x 250mm	110-150mm
<b>IN2-WR-C1032</b>	270mm	280mm	270mm x 250mm	90-170mm
<b>IN2-WR-B1033</b>	320mm	330mm	320mm x 250mm	140-250mm



For more information, call  
Customer Service on 1300 616 366  
or your local Product Specialist