

**Be Enlightened**



**MATRIX SURGICAL**  
affordable **healthcare** solutions

**A proudly Australian-owned company**

**t. 1300 616 366**

**e. [customerservice@matsurg.com.au](mailto:customerservice@matsurg.com.au)**

**w. [www.matrixsurgical.com.au](http://www.matrixsurgical.com.au)**



New and Unique from Matrix

## The Rutilight® The Light Beyond Your Hands



### Features:

- Rutilight is an LED ring light for illumination right at the point of operation
- Rutilight allows practical illumination of the surgical field, regardless of surrounding light, enabling clear access to difficult and badly lit areas
- Rutilight is an LED light ring which can be adapted to virtually any surgical instrument without affecting instrument function or hindering the surgeon's line of sight
- Rutilight can be fitted to any surgical instrument with a diameter between 5mm to 12mm, with a power source capable of accommodating 3 lights simultaneously, if required.



For more information, call  
Customer Service on 1300 616 366  
or your local Product Specialist

19 September 2018

### More Features:

- Similarly, the Rutilight does not hinder the balance or handling characteristics of your instruments due to super light weight construction.
- Rutilight is designed by surgeons for surgeons and offers colour temperature and rendition for organs, vessels and tissue with low power consumption LED Technology.
- Rutilight is ideal for open Surgeries including:

|   |   |
|---|---|
| <b>COLOPROCTOLOGY:</b> ○<br>Colon, rectal, anal and endoanal surgery.                         | <b>GENERAL AND DIGESTIVE SURGERY:</b> ○<br>Hepato-bilio-pancreatic and other surgery.       |
| <b>ENDOCRINOLOGY:</b> ○<br>Thyroid, parathyroid and suprarenal surgery and ganglion drainage. | <b>PAEDIATRIC SURGERY:</b> ○<br>Solid tumours, severe malformations and stenosis.           |
| <b>BREAST SURGERY:</b> ○<br>Subcutaneous areola-sparing mastectomy and underarm surgery.      | <b>PLASTIC AND RECONSTRUCTIVE SURGERY:</b> ○<br>Complex maxillofacial and other operations. |
| <b>OESOPHAGOGASTRIC SURGERY:</b> ○<br>Abdominal and cervical phases.                          | <b>TRAUMATOLOGY:</b> ○<br>Complex shoulder, hip, knee, foot and ankle injuries.             |
| <b>MORBID OBESITY:</b> ○<br>Open surgery.   | <b>UROLOGY:</b> ○<br>Pathologies affecting the urinary tract and reproductive system.       |
| <b>TRANSPLANTS:</b> ○<br>Liver, kidney, pancreas and lung.                                    | <b>GYNAECOLOGY:</b> ○<br>Hysterectomy, abdominal myomectomy and complex surgery.            |
| <b>PERITONEAL ONCOLOGY</b> ○<br>Peritoneal carcinomatosis.                                    |   |

- Rutilight is supplied as sterile single-use, in a box of 25 individual packages, and naturally, Rutilight is both CE marked and TGA approved.

TGA Certificate : 307598

### PRODUCT CODES

|                |   |
|----------------|---|
| <b>RL-01</b>   | Rutilight, sterile single-use box of 25 |
| <b>FA-RL01</b> | 3 port reusable power source            |

***The Rutilight provides from 25,000 lux at 5cm up to 90,000 lux at 2cm***

*"Enhancing surgeons visions of the surgical field."*

**Ask your Matrix Product Specialist for a demonstration  
or call Matrix direct  
on 1300 616 366**

**MATRIX SURGICAL**  
affordable **healthcare** solutions