

Ligation Clips & Appliers

Features

- High-quality titanium clips with chevron-shape for complete vessel overlapping
- MRI compatible and radio opaque
- Everest appliers with diamond dust jaws
- Cross-serration gives very secure closure, preventing slippage
- Precise smart-loading mechanism
- Adhesive backing on clip cartridge
- Colour-coded appliers for easy identification to the matching clip size

PRODUCT CODES: CLIPS

ARTG 235137

MDD-EVC010	Small (Red) 6 clips/cartridge 30 cartridges/box
MDD-EVC020	Medium (Blue) 6 clips/cartridge 30 cartridges/box
MDD-EVC030	Medium-Large (Green) 6 clips/cartridge 20 cartridges/box
MDD-EVC040	Large (Orange) 6 clips/cartridge 20 cartridges/box
MDD-EVC060	Micro (White) 6 clips/cartridge 30 cartridges/box
MDD-EVC020C	Medium 15x24 clips 360 clips per box
MDD-EVC010C	Small 15x24 clips 360 clips per box

PRODUCT CODES: CLIP APPLIERS

ARTG 234381

MDD-CEVC111	Small (Red) 20cm, each(CARDIAC)
MDD-CEVC121	Medium (Blue) 20cm, each(CARDIAC)
MDD-EVC111	Small (Red) 20cm, each
MDD-EVC113	Small (Red) 28cm, each
MDD-EVC121	Medium (Blue) 20cm, each
MDD-EVC123	Medium (Blue) 28cm, each
MDD-EVC131	Medium-Large (Green) 20cm, each
MDD-EVC133	Medium-Large (Green) 28cm, each
MDD-EVC141	Large (Orange) 20cm, each
MDD-EVC161	Micro (White) 20cm, each
MDD-EVC160	Micro (White) 15cm, each



For more information, call
Customer Service on 1300 616 366
or your local Product Specialist

