

MATRIX SURGICAL

affordable healthcare solutions

A proudly Australian-owned company

t. 1300 616 366

e. customerservice@matsurg.com.au

w. www.matrixsurgical.com.au

GLUBRAN®2 Surgical Glue

Summary:

- Synthetic modified co-monomer cyanoacrylate glue
- TGA certified for internal and endovascular use

Properties:

- Strong adhesive
- Haemostatic
- Sealant
- Bacteriostatic
- Latex-free

Features:

- Easy preparation (no thawing, mixing or waiting)
- Fast polymerisation (beginning in 1-2 seconds and completed in 60-90 seconds)
- Low polymerisation temperature (45°C maximum)
- The product is fully synthetic
- Non-toxic, non-mutagenic and bio-compatible
- Fully rebatable for private patients
- Efficacy in a wet environment
- Wide range of applications, including surgical, endoscopic and interventional radiology

Presentation:

- Individually packaged
- Box of 10 x 1ml vials



Sterilisation:

Packaging is gamma irradiated.

Storage:

- Refrigerate between 2°C and 8°C (do not freeze)
- 2 year shelf life from date of manufacture

ARTG Number: 158438

Rebate code: MG038

PRODUCT CODES

GEM-G-NB-2	Synthetic Surgical Glue (10 x 1ml vials)
-------------------	--



For more information, call
Customer Service on 1300 616 366
or your local Product Specialist